

Klorkleen Pack Formats

Klorkleen tablets are available in 1670mg tablet size. US EPA registered tablet size and formulation also available. The tablets are presented in a tub with a 2 year shelf life.

How to use Klorkleen

The correct strength of solution is produced by placing a tablet into the following volumes of water:

Use	1 tablet	Solution Strength
Floors,walls and areas of general disinfection	5 litres of water	200ppm available chlorine
Food handling/processing surfaces/equipment	10 litres of water	100ppm available chlorine
High Risk Areas (eg Operating Theatres)	1 litre of water	1000ppm available chlorine

Note: Responsibility for control and issue of guidelines for disinfection rests with the relevant medical authority and infection control personnel.

Medentech[®]
SAFER WATER, SAFER WORLD

Clonard Road, Wexford, Ireland
Tel: +353 53 9117900 Fax: +353 53 9141271

Email: info@medentech.com

www.medentech.com

Klorkleen[®] is a registered trademark.



2-29-031B

Klorkleen[®]

Effervescent Disinfectant & Cleaning Tablets

Multi-purpose disinfection & cleaning



US EPA registered tablet size & formulation available.

High level broad-spectrum cleaning & disinfection

Two essential jobs in one

- Powerful disinfection and thorough cleaning



What is Klorkleen?

Klorkleen combines a detergent and disinfectant in each tablet making it a multipurpose cleaning and disinfectant tablet. It does two essential jobs in one tablet- powerful disinfectant and thorough cleaning.

Klorkleen kills germs fast & washes thoroughly. Surfaces and utensils are left clean and safe.

What is Klorkleen used for?

- High, intermediate and low level environmental surface disinfection.
- Washing and cleaning surfaces.
- Klorkleen is specially formulated to give safe and easy-to-prepare cleaning and disinfecting solutions, with broad spectrum activity against micro-organisms including:
 - **Bacteria:** Escherichia spp, Campylobacter, C.difficile, E.coli O157, Enterococcus, Leptospira, Listeria, M.R.S.A., Pseudomonas spp, Salmonella, Staphylococcus, Streptococcus.
 - **Mycobacteria:** Tuberculosis.
 - **Fungi:** Aspergillus, Candida and others.
 - **Spores:** Bacillus.
 - **Viruses:** Human Immunodeficiency virus (HIV), Hepatitis B virus (HBV), Herpes, Poliomyelitis, SARS coronavirus, Bird Flu, H1N1 (Swine Flu)



Where is Klorkleen used?

Klorkleen can be used in a wide range of areas including

- High level Hospital & Clinic disinfection (eg Operating Theatres and general wards)
- Catering outlets such as Hotels, Restaurants, Bars, Canteens
- Klorkleen kills Salmonella, E Coli and Listeria
 - Germs lurking in the kitchen, on tables and on floors.
- Commercial and Industrial
- Klorkleen is ideal for cleaning/disinfecting surfaces in food and drink processing plants, factories, offices and institutions of all kinds.
- Sports Clubs, Gyms, Swimming Pools, Hairdressers, Beauty Salons
 - Washing the floor with Klorkleen will eliminate the germs that cause athletes foot.

Klorkleen complements Klorsept as part of a comprehensive infection control programme.

- Now you can disinfect at the same time as washing at no extra cost or effort.

Klorkleen Advantages

- Klorkleen kills germs fast and washes thoroughly.
- There is no need to measure, pour out or mix liquids, just drop a Klorkleen tablet in water and it dissolves fast to give the correct solutions.
- Tablets don't leak, spill or splash!
- Tablets are light and safe to handle.
- Klorkleen packs are compact and therefore use little storage space.
- With just one tub of Klorkleen 150 buckets of powerful cleaning and disinfectant solution is achieved.
- Will not stain or leave residue on high gloss floors
- Special low foam / high clean formulation
- US EPA registered tablet size & formulation available.